

## Sample field devices

### Spine Folders: PD series

The sample making push you to search for absolute quality, not only for one piece, but for the whole production.

In the meantime you are pressed to deliver quicker and quicker.

A lot of routine operations oblige you to take your time and to pay attention to keep the quality high.

These are often operations that simple tools could simplify.

The *Spine folder* of **PD** series try to answer to these needs.



PD1700

The *Spine folder* of **PD** series help you to make spine folder with hot melt glue.

They let you make an uniform pressure over the whole surface and lasting it the time need by the glue.




PD2700

The *Spine folder* of **PD** series have a completely manual and easy use.  
The **PD** series is supplied in two way:

- ◆ one place, **PD 1000** type. this solution is useful with double size Spine, up to 700mm
- ◆ two places, **PD 2000** type, shape. This solution is made for spine wide up to 320mm it is useful to prepare two spine simultaneously, so you can prepare the next spine while in the first spine the glue get cold.

In the picture you can see the standard supply and all its parts included.

## Technical Characteristics



Model PD1700



Model PD 2700

Description	PD 1700	PD 2700	PD 2920
Pressing elements	1	2	2
Laying surface			
depth		400mm	
width	700mm		920mm
Pressing Surface			
depth	40mm	40mm	
width	680mm	320mm	
Spine height	up to 40mm with 20 mm blocks up to 60mm with 40 mm blocks		

By request	PD 1700	PD 2700	PD 2920
Blocks for higher spines		By request	
Single pressing elements	-	yes	yes
Pressing Surface			
depth	60mm	60mm	60mm
width	680mm	320mm	860mm
Spine height	20 mm over block height		

**JEM** s.r.l.  
via F. Mola, 17  
20156 Milan  
Italy

Phone.: +39 02 48 49 80 16  
email: [info@jem2009.eu](mailto:info@jem2009.eu)  
site: [www.jem2009.eu](http://www.jem2009.eu)

---

